BAUSCH + LOMB

Oxane® 1300/5700

Silicone Oil

EN

Product Information

COMPOSITION 1 ml of highly pure sterile silicone oil, 1000-1500 mPa·s or 5000-5900 mPa·s.
Refractive index = 1.40.

INDICATIONS: Oxane® is used for intraocular tamponade in severe cases of retinal detachment, e.g. severe detachment with massive proliferative vitreo-retinopathy, traumatic detachment, giant tears as well as other retinal detachments, which cannot be treated with other forms of therapy. The use of silicone oils in cases of central foramina with detachment or in cases of severe diabetic retinopathies is at present under discussion. **CONTRAINDICATIONS:** Intraocular lenses made from silicone can be opacified by silicone oil. Therefore, silicone oil should not be applied in eyes with silicone IOL's. Other contraindications are not known. **METHOD OF USAGE AND DOSAGE:** Oxane® is injected intraocularly in a dose to be determined by the ophthalmic surgeon. The selection of patients, the choice of surgical procedure and the length of the indwelling period of the silicone oil in the eye are the sole responsibility of the ophthalmic surgeon. **PRECAUTIONS OF USE:** For ocular use only; do not sterilise; do not use if sterile barrier is breached; verify the expiry date; do not overfill the posterior segment with Oxane®; do not reuse; Oxane® is free from preservatives and can become contaminated if the single dose vial is used more than once or for multiple patients, any left-over fluid must be disposed. Infection risk is not just limited to bloodborne pathogens. Improper use of single-dose/single-use vials can result in lifethreatening bacterial infections including bloodstream infections, meningitis, and epidural abscesses. The stability and compatibility of silicone oils with ophthalmic laser and cryotherapy has not been studied. **WARNINGS:** Heavy liquids like all types of fluorocarbon liquids have to be removed completely from the vitreous cavity before silicone oil can be filled into the eye. Following the use of silicone oils that were not purified, the occurrence of corneal decompensation, cataract-highly act, secondary glaucoma and emulsification has been reported. Although the casual relationship with the purity of the silicone oil has not yet been established, the incidence of these complications appears to be lower following the use of highly pure silicone oil. The duration of the indwelling period of the silicone oil in the eye must be determined by the treating physician. Some reports indicate a period of use of more than two years. It has been shown that on removal of the silicone oil in the post-operative phase, during which the visual acuity neither increases nor further decreases again, the final operative result can be improved.

INTERACTION WITH OTHER SUBSTANCES: Pre interactions between Oxane® and other substances **STORAGE INSTRUCTIONS:** Store Silicone Oil below 25°C. **PRESERVATION AND STERILISATION:** Oxane® is supplied in a 10 mL glass syringe. Sterilisation: product and outer packaging simultaneously in a steam steriliser. **ADVERSE REACTIONS:** Adverse reactions and / or potentially sight-threatening complications that may be reasonably regarded as related and that were not previously expected in nature or degree of incidence should be reported to a Bausch & Lomb representative.

MEMBER
GREEN DOT
SCHEME

STERILE
STERILIZED USING STEAM

UPPER
TEMPERATURE
LIMIT

Bausch & Lomb Incorporated,
1400 North Goodman Street,
Rochester, NY 14609, USA

EC REP

Bausch & Lomb GmbH,
Brunnstrasse 1 Dam 165-173,
13581 Berlin, Germany

®/TM denotes trademarks of Bausch & Lomb Incorporated or its affiliates.

© 2015 Bausch & Lomb Incorporated. All rights reserved.

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®/TM denotes trademarks of Bausch & Lomb or its affiliates.
®/TM denotes trademarks of Bausch & Lomb Inc. or its affiliates.
©/TM denotes trademarks of Bausch & Lomb and/or its affiliates.