ENGLISH

INSTRUCTIONS FOR USE

Follow instructions for use, manufacturer guidelines, and institution procedures for flush administration.

DESCRIPTION
BD PosiFlush™ XS Syringe is a ready to use sterile medical device (according to regulation (UE) 2017/745 of the European Parliament and of the Council). It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. The contents of our unopened or undamaged blister packages are guaranteed to be sterile, non-toxic, and non-pyrogenic.

INTENDED USE/INDICATIONS FOR USE

- BD PosiFlush™ XS Syringes are intended to be used FOR FLUSHING ONLY of in-situ peripheral intravenous catheters (PIVC), peripherally inserted central catheters (PICCs), central venous catheters (CVCs), and implanted venous access ports.
- BD PosiFlush™ XS Syringe is not intended for dry product reconstitution, for medication dilution, or where intravenous therapy with sodium chloride is indicated.

- Using aseptic technique, **BD PosiFlush™ XS Syringe can be used on a sterile field.**

INTENDED PATIENT POPULATION

BD PosiFlush™ XS Syringe is to be used with patients with in-situ peripheral intravenous catheters (PIVCs), peripherally inserted central catheters (PICCs), central venous catheters (CVCs), and implanted venous access ports.

INTENDED USER

BD PosiFlush™ XS Syringe is to be used by healthcare professionals experienced in vascular access and the use of these devices.

The BD PosiFlush™ XS Syringes are manual devices that may be operated by a large range of people with various human characteristics, including hand sizes and strengths.

CONTRAINDICATIONS

- Use in hypernatremic and fluid retention patients, when the administration of sodium or chloride could be clinically detrimental.
- Do not use BD PosiFlush™ XS Syringe if a patient has a known allergy to any of its components, materials or 0.9% sodium chloride solution, which may lead to an allergic response in anaphylaxis.

PERFORMANCE CHARACTERISTICS

- Interoperability with Vascular Access Devices (VAD).
- All syringes (3 mL, 5 mL, and 10 mL) have the same 10 mL diameter syringe barrel and therefore the flush pressure is equivalent for all sizes.

WARNINGS

- Do not use if unit package or content is damaged.
- Do not use if product has been left at freezing temperature.
- Verify the expiration date on the product package or label. Do not use if product has expired.
- Do not use if syringe tip cap or stopper is damaged.
- Do not use if solution is cloudy or colored, contains a precipitate, or has any type of suspended particulate matter.
- Do not re-use. Re-use may lead to infection or other illness/injury.
- Small parts are a potential choking hazard. After use, discard small parts according to your facility protocol.
- Do not re-sterilize.
- Possible complications and/or adverse reactions associated with flushing may include sepsis, infections (local/systemic), mucocutaneous blood exposure, exposure to bloodborne pathogens, air embolism, particulate embolism, blood clots, phlebitis, leakage that may lead to hazardous drug/fluid exposure, irritation, or transitory taste or odor.
- Single blister pack consists of paper and Polypropylene - Polyamide film.

CLINICAL BENEFITS

- The BD PosiFlush™ XS Syringe is a pre-filled, single use 0.9% sodium chloride syringe that helps to improve clinical efficiency by eliminating steps and time involved in the manual preparation of saline syringes.
- The BD PosiFlush™ XS syringe reduces risk of touch contamination that may occur during manual preparation of saline flush syringes.

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