## ENGLISH

### INSTRUCTIONS FOR USE

**1. DESCRIPTION**  
BD PosiFlush™ SP Syringe is a ready to use medical device (according to 93/42/EEC medical device directive). It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. The syringe contents is guaranteed to be sterile, non-toxic and non-pyrogenic.

**2. INTENDED USE**  
- BD PosiFlush™ SP Syringe is intended **FOR FLUSHING ONLY** in-situ vascular access devices.

- BD PosiFlush™ SP Syringe is not intended for dry product reconstitution, for medication dilution, or where intravenous therapy with sodium chloride is indicated.

- **BD PosiFlush™ SP Syringe must not be used on a sterile field.**

### 3. CONTRAINDICATIONS

### 4. CAUTIONS

Do not re-sterilise before use.  
Do not use if unit package or content is damaged.  
Verify the expiration date on the product package or label. Do not use if product has expired.  
Do not use if syringe tip cap or stopper is damaged in any way that suggests or indicates syringe leakage.

Do not use if solution contains a precipitate, or has any type of suspended particulate matter. For single use only. Discard any partially used product.

Store at controlled temperature (15-25°C). Do not leave at freezing temperature.

Federal (USA) Law restricts this device to sale by or on the order of a physician.

Re-use may lead to infection or other illness/injury.

**5. DIRECTIONS FOR USE**  
Single use, single patient device.  
Use aseptic technique throughout the procedure.

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5. Push syringe plunger to expel the air. (Fig. 3)

6. Connect BD Luer-Lok™ syringe to vascular access device, taking care that there is no touch contamination of the connection. (Fig. 4)

7. Push syringe plunger to flush the required volume of saline following institution's policy. (Fig. 5)

In case of extreme plunger resistance, it is recommended that excessive force is not exerted. Some patients may experience a transitory taste or odour during flushing. This minor effect ceases shortly after the procedure is completed.

8. After use, dispose of in accordance with recognised procedure in your institution. (Fig. 6) DO NOT RE-USE.

**6. COMPOSITION PER UNIT**  
Polypropylene syringe with BD Luer-Lok™ tip. Elastomeric stopper free of natural rubber latex. Saline solution: sodium chloride 9 g/l (NaCl 0.9%), distilled sterile water to volumes, preservative free and non-pyrogenic.

Flow pack consists of Polypropylene.

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BD, the BD Logo and PosiFlush are trademarks of Becton, Dickinson and Company. © 2018 BD and its subsidiaries. All rights reserved.